

Riviera / Lounge-chair
Design / Lucidi / Pevere

Test
EN581-1 / EN581-2

Designed by Lucidi/Pevere, Riviera is a collection consisting of a chair, small armchair, lounger and stool; its main characteristic is an elegant minimal graphic design. The entirely welded structure is made from Fe360 tubular steel and rounded (half almond) bars. EMU-Coat finish (cataphoresis + painting).




Description

Code: **437**
Typology: **Armchairs**
Collection: **Riviera**

Stackable: 8
Static Load: 200 Kg
Cover: Yes
Cover Type: FC/1216
Test: EN581-1 / EN581-2
Cushion Material: Tessuto Standard

Technical Info



Width: **66.5 cm**
Depth: **72 cm**
Height: **69 cm**
Weight: **9.20 Kg**
Static Load: **200 Kg**

Packaging

Packages: 1
Pieces for pack: 2
Dimensions: **82x75x69 cm**
Volume: **0.424 mc**

Polyethylene bag and cardboard box

Colors

16
Marine Blue23
Matt White33
Pink Magnolia43
Sapphire Blue75
Dark Green17
Military Green24
Black37
Grey/Green46
Intense Red22
Antique Iron26
Maple Red41
Indian Brown62
Curry Yellow

Riviera / Lounge-chair
Design / Lucidi / Pevere

Test
EN581-1 / EN581-2

Materials

Steel

Iron and carbon alloy, with carbon content below 2%, treated to withstand the elements with the unique EMU-Coat anti-corrosio process.

Maintenance: To keep the product in good condition for a long time, we recommend storing it indoors and in a dry place during the winter to prevent condensate from forming. Before the winter and on a quarterly basis, if the products are kept near the sea, it is recommended to clean the metal surfaces with a soft cloth using water or detergent and protect them with Vaseline oil or car wax.

Cushion Material

Tessuto Standard

Maintenance: To keep the product in good condition for a long time, we recommend storing it indoors and in a dry place during the winter to prevent condensate from forming. Before the winter and on a quarterly basis, if the products are kept near the sea, it is recommended to clean the metal surfaces with a soft cloth using water or detergent and protect them with Vaseline oil or car wax.

Cushions for this collection: C/437 Seat cushion